


PIRONHA

FOOD IMPACTION TREATMENT SYSTEM



CORE
through the
impaction

REMOVE
the food

CLEAR
the field

PiranhaMedical.com

The **Piranha GI** catheter is advanced through an endoscope



DESIGN AND FUNCTIONALITY

Proximal 5" of shaft to add stiffness / pushability

Custom molded handle to interface with vacuum

Molded grips for ease of handing and manipulation

Internal liner is Teflon coated providing a lubricious surface to facilitate food suction

Custom tip with beveled edges and increased stiffness aides in coring through the impaction